



4.46 - High Resolution Laser spectroscopy

▷ Diode laser characteristics:




simple diode laser: p-n junction in forward bias: (homojunction)



better: double heterojunction → smaller active region → high divergence and also astigmatism

LNTR configuration → smaller linewidth




Tuning of output frequency:

- Laser temperature (slow)
- Laser current (fast)
- grating angle (only in single mode)

$$\Delta \omega^2 + \Delta \omega^2 \cdot \frac{\kappa}{P_{sat}} = 0$$

$$\Rightarrow \kappa = -P_{sat}$$

▷ Laser power and intensity:



$$\Gamma(L) = e^{-\chi \frac{1}{1+s} \cdot L} \approx 1 - \alpha \frac{1}{1+s}$$

$$R_{ing} \propto R = 1 - \Gamma$$

$$= \frac{\alpha}{1+s}$$

▷ Fabry - Perot interferometer

periodic transmission: $\Gamma_{FPI} = \frac{1}{1 + F \sin^2 \frac{\phi}{2}}$


free spectral range: $\Delta \omega_{FSR} = \frac{c}{2nL}$

$$\frac{\omega}{c} = n k_0 L = \frac{n \omega L}{c}$$

$$F = \frac{4R}{(1-R)^2}$$

FWHM = $\frac{\Delta \omega_{FSR}}{F}$ $\Gamma = \frac{\pi R}{1-R}$ considering also transv. modes: $\Delta \omega = \frac{c}{2nL} = \frac{c}{4nL}$

▷ Properties of rubidium Notation: $n^{2s+1}L_J$



▷ linear spectroscopy

- ~~LABERT-Beer~~: $\Gamma(\omega) = e^{-\chi(\omega)L}$ $\chi(\omega) = \chi_0 \frac{(\Delta\omega/2)^2}{(\omega - \omega_0)^2 + (\Delta\omega/2)^2}$ (LORENTZ)

- ~~DOPPLER~~-width: $\chi_0(\omega) = \chi_0(\omega_0) \exp\left[-\frac{(\omega - \omega_0)^2}{2\Delta\omega_D^2}\right]$

↳ obtained by convoluting $\chi(\omega)$ with velocity distribution

↳ natural width
 $\tau = \frac{1}{\Gamma} = \frac{1}{\Delta\nu_{\text{DOPPLER}}}$

▷ non-linear spectroscopy

- saturation spectroscopy: significant fraction of atoms is in excited state

$\chi(\omega) \rightarrow \tilde{\chi}(\tilde{\omega})$ $\Delta\tilde{\omega} = \omega \sqrt{1+s}$

absorbed power: $\Delta P \propto \frac{s}{1+s} \rightarrow \infty$ for $s \rightarrow \infty$ "optical power broadening"

pump beam: turn ~~LABERT~~-like into absorption spectrum

probe beam: measure the absorption on resonance \rightarrow ~~LORENTZ~~-Dip

(cross-over resonances \rightarrow as name suggests)

$E = h \cdot \omega = h \cdot \frac{c}{\lambda}$

$\omega = \frac{c}{\lambda}$

$\lambda = \frac{c}{\omega}$

